

Quality assurance monitoring, the lynchpin of rapid HIV screening in 10 publicly funded New York City sexually transmitted disease clinics, NYC, January, 2004

Kowalski, AV, Mavinkurve, M, Paneth-Pollak, R, Borrelli, J, Rubin, S, Blank, S Bureau of Sexually Transmitted Disease Control, New York City Department of Health and Mental Hygiene

PROJECT

- The New York City Department of Health and Mental Hygiene's Bureau of Sexually Transmitted Disease Control (NYCDOHMH BSTDC) operates 10 STD clinics throughout New York City.
 - Services are free and offered to individuals 12 vears and older, without parental consent, proof of insurance or citizenship.
 - Irrespective of chief complaint, all patients are offered HIV testing and counseling (confidential or anonymous).
- Approximately 500 cases of HIV are diagnosed annually in BSTDC STD clinics.
- In January of 2004, BSTDC implemented the OraQuick Advance Rapid HIV-1/2 Antibody test
 - HIV test results are provided within 20 minutes of specimen collection.
- OraQuick Advance Rapid HIV-1/2 Antibody test provides either a final negative result or a preliminary positive result (confirmed via Western Blot testing at the NYC Public Health Laboratory).

ISSUE

- Operational issues were addressed prior to implementing OraQuick Advance Rapid HIV-1/2 Antibody testing in the clinics
 - 1) Developed Rapid HIV testing protocols for staff and laboratory permits
 - 2) Re-engineered clinic workflows to accommodate high volume testing
 - 3) Trained and re-assigned staff to conduct testing, as conventional EIA screening had been previously performed in a reference laboratory
 - 4) Customized patient materials and pre/post test counseling messages for Rapid HIV testing
 - 5) Instructed staff to monitor and report unexpected results (e.g. preliminary positive, followed by western blot negative)
- A monthly average of 5 false-positives per 3900 tests was established as the false-positive baseline.
- A false-positive cluster was identified by BSTDC STD clinic staff in November, 2005.

RESULTS

40

- Prior to the introduction of OraQuick Advance Rapid
 - During the roll-out of rapid testing, 38,297 HIV tests were conducted, of these 15.047 were rapid tests.
 - OraQuick Advance Rapid HIV-1/2 Antibody test

Rapid testing continued to increase HIV testing volumes in BSTDC STD clinics.

- During November, 2005 a cluster of false-positive (preliminary-positive OraQuick, Western Blot negative) tests was detected.
- · Quality assurance monitoring data gathered over the next two months showed a decline in test specificity.
- 99.07%: below the lower limit of FDA and manufacturer's specificity level of 99.8 (95% CI: 99.6%-99.9%)).
- A monthly average of 5 false-positives per 3900 tests was observed as the false-positive baseline.
- BSTDC took immediate actions to identify the source and modify practices in BSTDC STD clinics.
- · Investigations showed no correlation between falsepositive tests and kit lots, storage temperatures, or test processing.
 - BSTDC implemented immediate repeat testing using whole-blood (finger-stick) on all preliminary-positive OraQuick specimens.
 - Immediate counseling messages were adjusted based on the second test result.
- By December 10th 2005, the number of false-positive tests returned to acceptable levels (81/47.204 rapid tests: test specificity 99.82%).

LESSONS LEARNED

- 1. Aggressive quality assurance monitoring and good communication enabled BSTDC to preserve the capacity to offer rapid oral-fluid testina.
- 2. BSTDC adapted policy and protocols for rapid HIV testing.
- · BSTDC STD clinics follow all preliminarypositive cases with a whole-blood (finger-stick) rapid HIV test.
- Regardless, of the results of the second rapid HIV test (whole-blood), Western Blot testing is still conducted
- 2. BSTDC's ability to offer additional rapid wholeblood (finger-stick) HIV testing restored confidence in rapid oral-fluid test technology for staff and patients.
- 3. Rapid oral-fluid HIV testing continues to increase HIV testing volumes in BSTDC STD clinics.
- 4. BSTDC STD continues to monitor oral-fluid rapid testina.
 - Monthly data runs are created each month to highlight discordant HIV test results.
- Clinic staff is required to report discordances in HIV test results in BSTDC STD clinics.

ACKNOWLEDGMENTS

The authors would particularly like to thank those BSTDC STD employees without whose hard work and dedication this work could not have been possible.

2000 1500 1000 were being they the try the cast of they be Months **Discordant Results compared to Rapid HIV Tests**, 2005 - 2006

Rapid and Conventional HIV Tests,

NYC STD Clinics, 2004

and and a

Month

Number of False Positive Tests compared to

Rapid HIV Tests, 2005

ŝ

3000 Issue

₽ 2500

Conventional

5000 -

4500

4000

3500 .8 3000

≩ 2500

Rapid



Rapid HIV Tests - Discordant Results

- HIV-1/2 Antibody test (January, 2004), BSTDC performed approximately 33,000 HIV tests annually.
- In January of 2004, BSTDC implemented the

35 false-positive tests/3.753 tests (test specificity)